



DENTIRO® Light

Ready-to-use surface disinfectant for the rapid and residue-free disinfection of non-invasive medical devices. Perfume-free.



Intended Purpose



DENTIRO® Light is a perfume-free surface disinfectant for the rapid and residue-free disinfection and cleaning of non-invasive medical devices such as operating tables, gurneys, IV poles, and dental chairs. The low alcohol content of less than 50% ensures high material compatibility and reduces the potential for allergies and skin irritation. DENTIRO® Light is easily biodegradable and free of aldehydes.

Intended User



The product is intended for use by dental or healthcare professionals with specific knowledge about the disinfection of medical devices. The product is not to be used directly on patients.

Instructions for Use



For your protection, wear disposable gloves and safety goggles. Completely moisten and wipe the entire surface of the medical device. Leave disinfectant to act and adhere to contact time. Use on alcohol compatible materials only. Before first time use, test material compatibility by applying the disinfectant to a small inconspicuous area. A list of compatible and incompatible materials is available on our website. If you are uncertain as to the material composition of the medical device, contact the manufacturer before use. All serious accidents occurring in relation to the Medical Device should be reported to the manufacturer or the competent authority.

Limitation



FOR PROFESSIONAL USE ONLY. Not for use on invasive medical devices. Do not mix with other products. Make sure that all visible soiling has been removed before disinfection.

Handling



WARNING H226 Flammable liquid and vapour. H319 Causes serious eye irritation. P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P233 Keep container tightly closed. P280 Wear protective gloves/protective clothing/eye protection/face protection. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P337 + P313 If eye irritation persists: Get medical advice/attention. P501 Dispose of contents/container in accordance with local and national regulations. For more information, please refer to Safety Data Sheet available on our website.



Storage



Store in cool, dry and well-ventilated facility. Keep out from direct sunlight, heat, and open flame. Storage temperature: < 30°C.

Physical Properties



Physical state: Clear, non-viscous liquid	pH-value: Neutral
Colour: Colourless	Shelf life: 36 months (unopened container)
Scent: Of alcohol	

Composition in 100 g



100 g contains 31.3 g ethanol, 16.7 g isopropanol, auxiliaries and water.

Effective Spectrum



Bactericidal	EN 13727:2012+A2:2015*, EN 16615:2015^	60 sec
Tuberculocidal	EN 14348:2005*	60 sec
Mycobactericidal	EN 14348:2005^	60 sec
Yeasticidal	EN 13624:2013*, EN 16615:2015^	60 sec
Enveloped viruses	EN 14476:2013+A1:2015*	15 sec
Virucidal	EN 14476:2013 ^	60 sec

^ Tested under low organic load / * Tested under high organic load

Compatibility



Acetal polymer	Compatible
Acrylic glass	Compatible
Acrylonitrile butadiene styrene	Compatible
Aluminium	Compatible
Artificial leather	Compatible
Chlorinated polyvinyl chloride	Compatible
Polyamide	Compatible
Polybutylene terephthalate	Compatible
Polycarbonate	Compatible
Polyetherimide	Compatible
Polyethylene	Compatible
Polyethylene terephthalate	Compatible
Polyphenylene oxides	Compatible
Polyphenylsulfone	Compatible
Polypropylene	Compatible



Polystyrene	Compatible
Polysulfone	Compatible
Polyvinyl chloride	Compatible
Polyvinylidene fluoride	Compatible
Rhodium	Compatible
Silicone	Compatible
Thermoplastic elastomer	Compatible
Thermoplastic polyurethane	Compatible
Ultra-high-molecular-weight polyethylene	Compatible




Delivery Forms



CH-031005	250 ml bottle
CH-031013	1 litre bottle
CH-031025	5 litre canister

Legal Manufacturer, EU Representative, Importer, Distributor



Oro Clean Chemie AG • Allmendstrasse 21 • 8320 Fehraltorf • Switzerland
 +41(0)44 226 44 44 •  info@oroclean.com •  www.oroclean.com



Oro Clean Chemie s.r.o • Vinohradská 2828/151 • Žižkov • 130 00 Praha 3 • Czech Republic
 info@oroclean.cz •  www.oroclean.cz