



Surface disinfectant with 50% alcohol for the disinfection of medical devices.

actisolve® alcosept is an aldehyde-free, ready-to-use surface disinfectant based on ethanol and quaternary ammonium compounds for the rapid disinfection of surfaces of non-invasive medical devices such as operating tables.

INSTRUCTIONS FOR USE

For your protection, wear disposable gloves and safety goggles. Completely moisten and wipe the entire surface of the medical device. Leave disinfectant to act and adhere to contact time. Wipe medical devices, which come into contact with skin with a water-soaked single-use cloth to prevent the build-up of biocides. Use on alcohol compatible materials only. Before first time use, test material compatibility by applying the disinfectant to a small inconspicuous area. A list of compatible and incompatible materials is available on our website. If you are uncertain as to the material composition of the medical device, contact the manufacturer before use. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

COMPATIBILITY

| | |
|--|----------------|
| Acetal polymer | Compatible |
| Acrylic glass | Not compatible |
| Acrylonitrile butadiene styrene | Compatible |
| Artificial leather | Compatible |
| Chlorinated polyvinyl chloride | Compatible |
| Latex | Compatible |
| Polyamide | Compatible |
| Polybutylene terephthalate | Compatible |
| Polycarbonate | Compatible |
| Polyetherimide | Compatible |
| Polyethylene | Compatible |
| Polyethylene terephthalate | Compatible |
| Polyphenylene oxides | Compatible |
| Polyphenylsulfone | Compatible |
| Polypropylene | Compatible |
| Polystyrene | Compatible |
| Polysulfone | Compatible |
| Polyvinyl chloride | Compatible |
| Polyvinylidene fluoride | Compatible |
| Silicone | Compatible |
| Thermoplastic elastomer | Compatible |
| Thermoplastic polyurethane | Compatible |
| Ultra-high-molecular-weight polyethylene | Compatible |

HANDLING

DANGER H226 Flammable liquid and vapour. H318 Causes serious eye damage. H336 May cause drowsiness or dizziness. P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P233 Keep container tightly closed. P261 Avoid breathing fume/mist/vapours/spray. P280 Wear protective gloves/protective clothing/eye protection/face protection. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P501 Dispose of contents/container in accordance with local and national regulations. For more information, please refer to Safety Data Sheet available on our website.

EFFECTIVE SPECTRUM / EXPOSURE TIMES

| | |
|-------------------------------------|--------|
| Bactericidal | 60 sec |
| Antibiotic-resistant bacteria | 60 sec |
| Tuberculocidal | 60 sec |
| Mycobactericidal | 60 sec |
| Yeasticidal | 60 sec |
| Drug-resistant yeasts | 60 sec |
| Fungicidal | 2 min |
| Virucidal against enveloped viruses | 15 sec |
| Bovine viral diarrhea virus (BVDV) | 15 sec |
| Coronavirus | 15 sec |
| Hepatitis B virus | 15 sec |
| Hepatitis C virus | 15 sec |
| Herpes simplex virus | 15 sec |
| Human immunodeficiency virus (HIV) | 15 sec |
| Influenza A virus | 15 sec |
| Vaccinia virus | 15 sec |
| Virucidal | 60 sec |
| Adenovirus | 60 sec |
| Norovirus | 30 sec |
| Poliovirus | 60 sec |
| Polyomavirus | 30 sec |

Complies with EN 13624, EN 13727, EN 14476 and EN 16615.

COMPOSITION

100 g contains 50 g 1-propanol, 0.057 g N-alkyl aminopropyl glycine, auxiliaries and water.

PHYSICAL PROPERTIES

| | |
|-----------------|---------------------------|
| Physical state: | Clear, non-viscous liquid |
| Colour: | Colourless |
| Scent: | Of alcohol |
| pH-value: | Neutral |
| Shelf life: | 48 months |

DELIVERY FORMS

| | |
|-----------|------------------|
| AC-021005 | 250 ml bottle |
| AC-021060 | 1 litre bottle |
| AC-021065 | 5 litre canister |

STORAGE

Store in a cool, dry and well-ventilated place. Keep away from direct sunlight, heat, and open flame. Storage temperature: < 30 °C. Dispose of the product after expiry date.

LEGAL MANUFACTURER



United Disinfectant Manufacturers AG
Dr. Grass-Strasse 12
9490 Vaduz
Principality of Liechtenstein
Telephone: +423 237 15 03
info@udm.li
www.udm.li