



Disinfectant wipes for the residue-free disinfection of medical devices.

actisolve® alcowipes are ready-to-use disinfectant wipes for the residue-free disinfection of alcohol-resistant surfaces of non-invasive medical devices such as operating tables, treatment chairs, and medical equipment.

INSTRUCTIONS FOR USE

Write the use-by date after opening on the dispenser tub. For your protection, wear disposable gloves and safety goggles. Insert the refill pouch into the designated dispenser tub. Before first use, cut open the refill pouch. Pull the first wipe from the centre of the roll through the star-shaped opening in the lid. Close the dispenser tub firmly after each use. Completely moisten and wipe the entire surface of the medical device. Leave the disinfectant to act and adhere to contact time. Use on alcohol-compatible materials only. Before first use, check material compatibility on a small, inconspicuous area. A list of compatible and incompatible materials is available on our website. If you are unsure about the material composition of the medical device to be disinfected, contact the manufacturer before use. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

COMPATIBILITY

Acetal polymer	Compatible
Acrylic glass	Compatible
Acrylonitrile butadiene styrene	Compatible
Aluminium	Compatible
Artificial leather	Compatible
Chlorinated polyvinyl chloride	Compatible
Latex	Compatible
Polyamide	Compatible
Polybutylene terephthalate	Compatible
Polycarbonate	Compatible
Polyetherimide	Compatible
Polyethylene	Compatible
Polyethylene terephthalate	Compatible
Polyphenylene oxides	Compatible
Polyphenylsulfone	Compatible
Polypropylene	Compatible
Polystyrene	Compatible
Polysulfone	Compatible
Polyvinyl chloride	Compatible
Polyvinylidene fluoride	Compatible
Rhodium	Compatible
Silicone	Compatible
Thermoplastic elastomer	Compatible
Thermoplastic polyurethane	Compatible
Ultra-high-molecular-weight polyethylene	Compatible

EFFECTIVE SPECTRUM / EXPOSURE TIMES

Bactericidal	60 sec
Antibiotic-resistant bacteria	60 sec
Tuberculocidal	60 sec
Mycobactericidal	60 sec
Yeasticidal	60 sec
Drug-resistant yeasts	60 sec
Virucidal against enveloped viruses	15 sec
Bovine viral diarrhea virus (BVDV)	15 sec
Coronavirus	15 sec
Hepatitis B virus	15 sec
Hepatitis C virus	15 sec
Herpes simplex virus	15 sec
Human immunodeficiency virus (HIV)	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Virucidal	60 sec
Adenovirus	60 sec
Norovirus	60 sec
Poliovirus	60 sec
Rotavirus	30 sec

Complies with EN 13624, EN 13727, EN 14348, EN 14476 and EN 16615.

PHYSICAL PROPERTIES

Physical state:	Liquid absorbed on fleece
Colour:	Colourless
Scent:	Of alcohol
pH-value:	Neutral
Shelf life:	36 months

COMPOSITION

100 g impregnating solution contains 31.3 g ethanol, 16.7 g isopropanol, auxiliaries and water.

STORAGE

Store in a cool, dry and well-ventilated place. Keep away from direct sunlight, heat, and open flame. Storage temperature: < 30 °C. Dispose of the refill pouch 30 days after opening.



Disinfectant wipes for the residue-free disinfection of medical devices.

HANDLING

WARNING H226 Flammable liquid and vapour. H319 Causes serious eye irritation. P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P233 Keep container tightly closed. P280 Wear protective gloves/protective clothing/eye protection/face protection. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P337 + P313 If eye irritation persists: Get medical advice/attention. P501 Dispose of contents/container in accordance with local and national regulations. Further information can be found in the Safety Data Sheet available on udm.swiss.

DELIVERY FORMS

AC-041079	Refill pouch with 120 disinfectant wipes (14.5 x 20 cm), classic scent
AC-041084	Refill pouch with 120 disinfectant wipes (14.5 x 20 cm), lemon scent
AC-041092	Jumbo, Refill pouch with 100 disinfectant wipes (14.5 x 40 cm), classic scent
AC-041097	Jumbo, Refill pouch with 100 disinfectant wipes (14.5 x 40 cm), lemon scent
AC-042085	Saver, Refill pouch with 120 disinfectant wipes (14.5 x 20 cm), lemon scent
AC-042080	Saver, Refill pouch with 250 disinfectant wipes (14.5 x 20 cm), lemon scent
AC-042084	Soft, Refill pouch with 120 disinfectant wipes (19.3 x 21 cm), lemon scent

LEGAL MANUFACTURER



United Disinfectant Manufacturers AG
Dr. Grass-Strasse 12
9490 Vaduz
Principality of Liechtenstein
Telephone: +423 237 15 03
info@udm.li
www.udm.li

QMS-IFU-TDS-12

Rev.1.0 / 15.01.2026