




Concentrate for the disinfection, cleaning and deodorisation of dental suction units.

actisolve® aspisept is a concentrate for the disinfection and cleaning of dental suction units and spittoon bowls.

INSTRUCTIONS FOR USE

For your protection, wear a laboratory gown, household gloves and safety goggles. Disinfect the suction unit daily or in accordance with local regulations. Rinse the suction unit thoroughly with water after surgical procedures and prior to disinfection. Prepare 2 litres of working solution according to the dosing table. Aspirate 1.5 litres of the working solution and pour the remaining 500 ml down the spittoon bowl. Leave to act. Rinse thoroughly with water after disinfection. Not compatible with aluminium, copper, galvanized steel, and polyurethane. A list of compatible and incompatible materials is available on our website. If you are unsure about the material composition of the medical device to be disinfected, contact the manufacturer before use. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

DOSING

	Concentrate	+	Water	Contact Time
2%	40 ml		1960 ml	60 min

COMPATIBILITY

Acrylonitrile butadiene styrene	Compatible
Latex	Compatible
Polyethylene	Compatible
Polypropylene	Compatible
Polystyrene	Compatible
Polyvinyl chloride	Compatible

HANDLING

DANGER H290 May be corrosive to metals. H314 Causes severe skin burns and eye damage. H410 Very toxic to aquatic life with long lasting effects. P280 Wear protective gloves/protective clothing/eye protection/face protection. P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P391 Collect spillage. P501 Dispose of contents/container in accordance with local and national regulations. For more information, please refer to Safety Data Sheet available on our website.

COMPOSITION

100 g contains 2.475 g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine, 7.5 g quaternary ammonium compounds, anti-foaming agents, auxiliaries and water.

EFFECTIVE SPECTRUM / EXPOSURE TIMES

Bactericidal	2% in 5 min
Antibiotic-resistant bacteria	2% in 5 min
Yeasticidal	2% in 5 min
Drug-resistant yeasts	2% in 5 min
Virucidal against enveloped viruses	2% in 5 min
Bovine viral diarrhea virus (BVDV)	2% in 1 min
Coronavirus	2% in 5 min
Hepatitis B virus	2% in 5 min
Hepatitis C virus	2% in 1 min
Herpes simplex virus	2% in 5 min
Human immunodeficiency virus (HIV)	2% in 5 min
Influenza A virus	2% in 1 min
Vaccinia virus	2% in 5 min
Virucidal	2% in 60 min
Adenovirus	2% in 30 min
Norovirus	2% in 60 min
Poliovirus	2% in 30 min

Complies with EN 13624, EN 13727, EN 14476, EN 14561, EN 14562 and EN 17111.

PHYSICAL PROPERTIES

Physical state:	Clear, slightly viscous liquid
Colour:	Yellow
Scent:	Aromatic
pH-value:	10.0-12.0
Shelf life:	48 months

STORAGE

Store in a cool, dry and well-ventilated place. Keep away from direct sunlight, heat, and open flame. Storage temperature: < 30 °C. Dispose of the product after expiry date.

DELIVERY FORMS

AC-071020	2 litre bottle
AC-071025	5 litre canister

LEGAL MANUFACTURER



United Disinfectant Manufacturers AG
Dr. Grass-Strasse 12
9490 Vaduz
Principality of Liechtenstein
Telephone: +423 237 15 03
info@udm.li
www.udm.li