







Concentrate for the disinfection of dental impressions.

actisolve® impresoak is a concentrate for the disinfection of dental impressions made of alginate, polyether, polysulfide, A-silicone, and C-silicone. actisolve® impresoak does not affect the dimensional stability of dental impressions.

INSTRUCTIONS FOR USE

For your protection, wear disposable gloves and safety goggles. Prepare the working solution according to the dosing table. Immerse the dental impression and impression tray completely in the disinfectant solution. Adhere to, but do not exceed, the indicated contact time. After disinfection, rinse the dental impressions thoroughly under running water of at least drinking-water quality for 10 seconds. Replace the working solution daily or when visibly contaminated. A list of compatible and incompatible materials is available on our website. If you are unsure about the material composition of the medical device to be disinfected, contact the manufacturer before use. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

DOSING

	Concentrate 	+	Water 	Contact Time 
3%	30 ml		970 ml	3 min

COMPATIBILITY

Acrylic glass	Compatible
Acrylonitrile butadiene styrene	Compatible
Alginate	Compatible
A-silicone	Compatible
C-silicone	Compatible
Latex	Compatible
Polycarbonate	Compatible
Polyether	Compatible
Polyethylene	Compatible
Polyethylene terephthalate	Compatible
Polypropylene	Compatible
Polystyrene	Compatible
Polysulfide	Compatible
Polyvinyl chloride	Compatible
Stainless steel	Compatible
Thermoplastic elastomer	Compatible
Thermoplastic polyurethane	Compatible

HANDLING

DANGER H302 Harmful if swallowed. H314 Causes severe skin burns and eye damage. H410 Very toxic to aquatic life with long lasting effects. P280 Wear protective gloves/protective clothing/eye protection/face protection. P301 + P312 IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P501 Dispose of contents/container in accordance with local and national regulations. Further information can be found in the Safety Data Sheet available on udm.swiss.

EFFECTIVE SPECTRUM / EXPOSURE TIMES

Bactericidal	3% in 3 min
Antibiotic-resistant bacteria	3% in 3 min
Yeasticidal	3% in 3 min
Drug-resistant yeasts	3% in 3 min
Virucidal against enveloped viruses	3% in 3 min
Bovine viral diarrhea virus (BVDV)	3% in 3 min
Coronavirus	3% in 3 min
Hepatitis B virus	3% in 3 min
Hepatitis C virus	3% in 3 min
Herpes simplex virus	3% in 3 min
Human immunodeficiency virus (HIV)	3% in 3 min
Influenza A virus	3% in 3 min
Vaccinia virus	3% in 3 min

Complies with EN 13624, EN 13727, EN 14476, EN 14561, EN 14562 and EN 17111.

COMPOSITION

100 g contains 30 g benzalkonium chloride, auxiliaries and water.

PHYSICAL PROPERTIES

Physical state:	Clear, slightly viscous liquid
Colour:	Blue
Scent:	Aromatic
pH-value:	10.5-12.5
Shelf life:	36 months

STORAGE

Store in a cool, dry and well-ventilated place. Keep away from direct sunlight, heat, and open flame. Storage temperature: < 30 °C. Dispose of the product after expiry date.

DELIVERY FORMS

AC-071041	1 litre bottle
-----------	----------------

LEGAL MANUFACTURER



United Disinfectant Manufacturers AG
Dr. Grass-Strasse 12
9490 Vaduz
Principality of Liechtenstein
Telephone: +423 237 15 03
info@udm.li
www.udm.li