




Disinfectant concentrate for sensitive surfaces of medical devices.

actisolve® maxsurface is a disinfectant concentrate for the disinfection of large surfaces of non-invasive medical devices made of sensitive materials, such as radiology equipment. It is compatible with acrylic glass, polycarbonate (PC), and PVC and is free of aldehydes, phenols, and caustic potash.

INSTRUCTIONS FOR USE

For your protection, wear a laboratory gown, household gloves and safety goggles. Prepare the working solution according to the dosing table. Completely wipe or fully immerse the surface of the medical device to be disinfected and leave to act. Adhere to the indicated contact time. If necessary, dry the surface with a clean single-use cloth. Before first use, check material compatibility on a small, inconspicuous area. A list of compatible and incompatible materials is available on our website. If you are unsure about the material composition of the medical device to be disinfected, contact the manufacturer before use. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

DOSING

	Concentrate	+	Water	Contact Time
2%	20 ml		980 ml	5 min

COMPATIBILITY

Acetal polymer	Compatible
Acrylic glass	Compatible
Acrylonitrile butadiene styrene	Compatible
Aluminium	Compatible with restriction
Artificial leather	Compatible
Chlorinated polyvinyl chloride	Compatible
Latex	Compatible
Polyamide	Compatible
Polybutylene terephthalate	Compatible
Polycarbonate	Compatible
Polyetherimide	Compatible
Polyethylene	Compatible
Polyethylene terephthalate	Compatible
Polyphenylene oxides	Compatible
Polyphenylsulfone	Compatible
Polypropylene	Compatible
Polystyrene	Compatible
Polysulfone	Compatible
Polyvinyl chloride	Compatible
Polyvinylidene fluoride	Compatible
Silicone	Compatible
Thermoplastic elastomer	Compatible
Thermoplastic polyurethane	Compatible
Ultra-high-molecular-weight polyethylene	Compatible

HANDLING

WARNING H315 Causes skin irritation. H410 Very toxic to aquatic life with long lasting effects. P280 Wear protective gloves/protective clothing/eye protection/face protection. P302 + P352 IF ON SKIN: Wash with plenty of soap and water. P332 + P313 If skin irritation occurs: Get medical advice/attention. P362 Take off contaminated clothing. P391 Collect spillage. P501 Dispose of contents/container in accordance with local and national regulations.

Further information can be found in the Safety Data Sheet available on udm.swiss.

EFFECTIVE SPECTRUM / EXPOSURE TIMES

Bactericidal	1% in 1 min
Antibiotic-resistant bacteria	1% in 1 min
Yeasticidal	1% in 1 min
Drug-resistant yeasts	1% in 1 min
Virucidal against enveloped viruses	1% in 5 min
Bovine viral diarrhea virus (BVDV)	1% in 5 min
Coronavirus	1% in 5 min
Hepatitis B virus	1% in 5 min
Hepatitis C virus	1% in 5 min
Herpes simplex virus	1% in 5 min
Human immunodeficiency virus (HIV)	1% in 5 min
Influenza A virus	1% in 5 min
Vaccinia virus	1% in 5 min
Limited spectrum virucidal	2% in 5 min
Adenovirus	2% in 1 min
Norovirus	2% in 5 min

Complies with EN 13624, EN 13727, EN 14476 and EN 16615.

COMPOSITION

100 g contains 10 g quaternary ammonium compounds, auxiliaries and water.

PHYSICAL PROPERTIES

Physical state:	Clear, slightly viscous liquid
Colour:	Green
Scent:	Lemon
pH-value:	6.0-8.0
Shelf life:	36 months

STORAGE

Store in a cool, dry and well-ventilated place. Keep away from direct sunlight, heat, and open flame. Storage temperature: < 30 °C. Dispose of the product after expiry date.

DELIVERY FORMS

AC-010005	250 ml bottle
AC-010059	1 litre bottle
AC-010060	2 litre bottle
AC-010065	5 litre canister

LEGAL MANUFACTURER

UDM UNITED DISINFECTANT MANUFACTURERS

United Disinfectant Manufacturers AG
Dr. Grass-Strasse 12
9490 Vaduz
Principality of Liechtenstein
Telephone: +423 237 15 03
info@udm.li
www.udm.li

QMS-IFU-TDS-18

Rev.1.0 / 15.01.2026