



## Alcohol-free surface disinfectant for medical devices.

actisolve® alcofree is a ready-to-use solution for the disinfection of surfaces of non-invasive medical devices such as treatment chairs with high-gloss surfaces that are not compatible with alcohol. The disinfectant is suitable for the disinfection of surfaces of medical devices made of plastic, acrylic glass, artificial leather and other alcohol-sensitive materials.

### INSTRUCTIONS FOR USE

For your protection, wear disposable gloves and safety goggles. Completely moisten and wipe the entire surface of the medical device. Leave the disinfectant to act and adhere to contact time. Wipe medical devices that come into contact with skin using a water-soaked single-use cloth to prevent the build-up of biocides. Before first use, check material compatibility on a small, inconspicuous area. A list of compatible and incompatible materials is available on our website. If you are unsure about the material composition of the medical device to be disinfected, contact the manufacturer before use. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### COMPATIBILITY

Acetal polymer	Compatible
Acrylic glass	Compatible
Acrylonitrile butadiene styrene	Compatible
Artificial leather	Compatible
Ceramic	Compatible
Chlorinated polyvinyl chloride	Compatible
Latex	Compatible
Polyamide	Compatible
Polybutylene terephthalate	Compatible
Polycarbonate	Compatible
Polyetherimide	Compatible
Polyethylene	Compatible
Polyethylene terephthalate	Compatible
Polyphenylene oxides	Compatible
Polyphenylsulfone	Compatible
Polypropylene	Compatible
Polystyrene	Compatible
Polysulfone	Compatible
Polyvinyl chloride	Compatible
Polyvinylidene fluoride	Compatible
Rhodium	Compatible
Silicone	Compatible
Stainless steel	Compatible
Steel	Compatible
Thermoplastic elastomer	Compatible
Thermoplastic polyurethane	Compatible
Ultra-high-molecular-weight polyethylene	Compatible

### HANDLING

No signal word according to Regulation (EC) 1272/2008. H412 Harmful to aquatic life with long lasting effects. P273 Avoid release to the environment. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P501 Dispose of contents/container in accordance with local and national regulations. Further information can be found in the Safety Data Sheet available on [udm.swiss](http://udm.swiss).

### EFFECTIVE SPECTRUM / EXPOSURE TIMES

Bactericidal	60 sec
Antibiotic-resistant bacteria	60 sec
Yeasticidal	60 sec
Drug-resistant yeasts	60 sec
Virucidal against enveloped viruses	15 sec
Bovine viral diarrhea virus (BVDV)	15 sec
Coronavirus	15 sec
Hepatitis B virus	15 sec
Hepatitis C virus	15 sec
Herpes simplex virus	15 sec
Human immunodeficiency virus (HIV)	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Rotavirus	30 sec

Complies with EN 13624, EN 13727, EN 14476 and EN 16615.

### COMPOSITION

100 g contains 0.4 g benzalkonium chloride, 0.4 g N-alkylaminopropyl glycine, auxiliaries and water.

### PHYSICAL PROPERTIES

Physical state:	Clear, non-viscous liquid
Colour:	Colourless
Scent:	Apple
pH-value:	8.2-8.8
Shelf life:	48 months

### STORAGE

Store in a cool, dry and well-ventilated place. Keep away from direct sunlight, heat, and open flame. Storage temperature: < 30 °C. Dispose of the product after expiry date.

### DELIVERY FORMS

AC-011004	200 ml foamer bottle
AC-011005	250 ml bottle
AC-011100	1 litre bottle
AC-011102	2 litre bottle
AC-011104	5 litre canister

### LEGAL MANUFACTURER



United Disinfectant Manufacturers AG  
Dr. Grass-Strasse 12  
9490 Vaduz  
Principality of Liechtenstein  
Telephone: +423 237 15 03  
[info@udm.li](mailto:info@udm.li)  
[www.udm.li](http://www.udm.li)