



Surface disinfectant with a low alcohol content for medical devices.

actisolve® alcolite is a ready-to-use solution with a low alcohol content for the disinfection of surfaces of non-invasive medical devices. The product can be safely used for the disinfection of sensitive medical devices such as treatment chairs with high-gloss surfaces. The disinfectant is compatible with the most commonly used materials in medical devices as well as acrylic glass and artificial leather.

INSTRUCTIONS FOR USE

For your protection, wear disposable gloves and safety goggles. Completely moisten and wipe the entire surface of the medical device. Leave the disinfectant to act and adhere to contact time. Wipe medical devices that come into contact with skin using a water-soaked single-use cloth to prevent the build-up of biocides. Use on alcohol-compatible materials only. Before first use, check material compatibility on a small, inconspicuous area. A list of compatible and incompatible materials is available on our website. If you are unsure about the material composition of the medical device to be disinfected, contact the manufacturer before use. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

COMPATIBILITY

| | |
|--|------------|
| Acetal polymer | Compatible |
| Acrylic glass | Compatible |
| Acrylonitrile butadiene styrene | Compatible |
| Artificial leather | Compatible |
| Chlorinated polyvinyl chloride | Compatible |
| Latex | Compatible |
| Polyamide | Compatible |
| Polybutylene terephthalate | Compatible |
| Polycarbonate | Compatible |
| Polyetherimide | Compatible |
| Polyethylene | Compatible |
| Polyethylene terephthalate | Compatible |
| Polyphenylene oxides | Compatible |
| Polyphenylsulfone | Compatible |
| Polypropylene | Compatible |
| Polystyrene | Compatible |
| Polysulfone | Compatible |
| Polyvinyl chloride | Compatible |
| Polyvinylidene fluoride | Compatible |
| Silicone | Compatible |
| Thermoplastic elastomer | Compatible |
| Thermoplastic polyurethane | Compatible |
| Ultra-high-molecular-weight polyethylene | Compatible |

HANDLING

DANGER H225 Highly flammable liquid and vapour. H319 Causes serious eye irritation. P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P233 Keep container tightly closed. P280 Wear protective gloves/protective clothing/eye protection/face protection. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P337 + P313 If eye irritation persists: Get medical advice/attention. P501 Dispose of contents/container in accordance with local and national regulations. Further information can be found in the Safety Data Sheet available on udm.swiss.

EFFECTIVE SPECTRUM / EXPOSURE TIMES

| | |
|-------------------------------------|--------|
| Bactericidal | 60 sec |
| Antibiotic-resistant bacteria | 60 sec |
| Yeasticidal | 60 sec |
| Drug-resistant yeasts | 60 sec |
| Virucidal against enveloped viruses | 15 sec |
| Bovine viral diarrhea virus (BVDV) | 15 sec |
| Coronavirus | 15 sec |
| Hepatitis B virus | 15 sec |
| Hepatitis C virus | 15 sec |
| Herpes simplex virus | 15 sec |
| Human immunodeficiency virus (HIV) | 15 sec |
| Influenza A virus | 15 sec |
| Vaccinia virus | 15 sec |
| Limited spectrum virucidal | 5 min |
| Adenovirus | 5 min |
| Norovirus | 5 min |

Complies with EN 13624, EN 13727, EN 14476 and EN 16615.

COMPOSITION

100 g contains 7.5 g ethanol, 15 g isopropanol, 0.76 g quaternary ammonium compounds, auxiliaries and water.

PHYSICAL PROPERTIES

| | |
|-----------------|---------------------------|
| Physical state: | Clear, non-viscous liquid |
| Colour: | Colourless |
| Scent: | Of alcohol |
| pH-value: | Neutral |
| Shelf life: | 24 months |

DELIVERY FORMS

| | |
|-----------|------------------|
| AC-032005 | 250 ml bottle |
| AC-032013 | 1 litre bottle |
| AC-032020 | 2 litre bottle |
| AC-032025 | 5 litre canister |

STORAGE

Store in a cool, dry and well-ventilated place. Keep away from direct sunlight, heat, and open flame. Storage temperature: < 30 °C. Dispose of the product after expiry date.

LEGAL MANUFACTURER



United Disinfectant Manufacturers AG
Dr. Grass-Strasse 12
9490 Vaduz
Principality of Liechtenstein
Telephone: +423 237 15 03
info@udm.li
www.udm.li