



## Surface disinfectant with a high alcohol content for the residue-free disinfection of medical devices.

actisolve® alcopro is a ready-to-use solution with a high alcohol content for the rapid disinfection of alcohol-resistant surfaces of non-invasive medical devices such as operating tables, treatment chairs, and medical equipment with alcohol-compatible surfaces.

### INSTRUCTIONS FOR USE

For your protection, wear disposable gloves and safety goggles. Completely moisten and wipe the entire surface of the medical device. Leave the disinfectant to act and adhere to contact time. Use on alcohol-compatible materials only. Before first use, check material compatibility on a small, inconspicuous area. A list of compatible and incompatible materials is available on our website. If you are unsure about the material composition of the medical device to be disinfected, contact the manufacturer before use. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### COMPATIBILITY

Acetal polymer	Compatible
Acrylic glass	Not compatible
Acrylonitrile butadiene styrene	Compatible
Artificial leather	Compatible
Chlorinated polyvinyl chloride	Compatible
Ethylene vinyl acetate	Compatible
High impact polystyrene	Compatible
Latex	Compatible
Liquid crystal polymer	Compatible
Polyamide	Compatible
Polybutylene terephthalate	Compatible
Polycarbonate	Compatible
Polyether ether ketone	Compatible
Polyetherimide	Compatible
Polyethylene	Compatible
Polyethylene terephthalate	Compatible
Polyethylene terephthalate glycol	Compatible
Polyactic acid	Compatible
Polyphenylene oxides	Compatible
Polyphenylene sulfide	Compatible
Polyphenylsulfone	Compatible
Polypropylene	Compatible
Polystyrene	Compatible
Polysulfone	Compatible
Polyurethane	Compatible
Polyvinyl chloride	Compatible
Polyvinylidene fluoride	Compatible
Silicone	Compatible
Thermoplastic elastomer	Compatible
Thermoplastic polyurethane	Compatible
Ultra-high-molecular-weight polyethylene	Compatible

### HANDLING

DANGER H226 Flammable liquid and vapour. H318 Causes serious eye damage. H336 May cause drowsiness or dizziness. P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P233 Keep container tightly closed. P261 Avoid breathing fume/mist/vapours/spray. P280 Wear protective gloves/protective clothing/eye protection/face protection. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P501 Dispose of contents/

container in accordance with local and national regulations. Further information can be found in the Safety Data Sheet available on [udm.swiss](http://udm.swiss).

### EFFECTIVE SPECTRUM / EXPOSURE TIMES

Bactericidal	60 sec
Antibiotic-resistant bacteria	60 sec
Tuberculocidal	60 sec
Mycobactericidal	60 sec
Yeasticidal	60 sec
Drug-resistant yeasts	60 sec
Virucidal against enveloped viruses	15 sec
Bovine viral diarrhoea virus (BVDV)	15 sec
Coronavirus	15 sec
Hepatitis B virus	15 sec
Hepatitis C virus	15 sec
Herpes simplex virus	15 sec
Human immunodeficiency virus (HIV)	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Limited spectrum virucidal	30 sec
Adenovirus	30 sec
Norovirus	30 sec

Complies with EN 13624, EN 13727, EN 14476 and EN 16615.

### COMPOSITION

100 g contains 25 g ethanol, 35 g 1-propanol, auxiliaries and water.

### PHYSICAL PROPERTIES

Physical state:	Clear, non-viscous liquid
Colour:	Colourless
Scent:	Of alcohol
pH-value:	Neutral
Shelf life:	24 months

### STORAGE

Store in a cool, dry and well-ventilated place. Keep away from direct sunlight, heat, and open flame. Storage temperature: < 30 °C. Dispose of the product after expiry date.

### DELIVERY FORMS

AC-033005	250 ml bottle
AC-033013	1 litre bottle
AC-033020	2 litre bottle
AC-033025	5 litre canister

### LEGAL MANUFACTURER



United Disinfectant Manufacturers AG  
Dr. Grass-Strasse 12  
9490 Vaduz  
Principality of Liechtenstein  
Telephone: +423 237 15 03  
[info@udm.li](mailto:info@udm.li)  
[www.udm.li](http://www.udm.li)