## EN (English)

**INTENDED PURPOSE:** ULTRACLEAN-3 is a concentrate for the cleaning of surgical and dental instruments in an ultrasonic cleaning device.

**INTENDED USERS:** The product is intended for use by dental and healthcare professionals who have knowledge of the disinfection of medical devices. The product must not be used directly on patients.

**INSTRUCTIONS FOR USE:** For your protection, wear a laboratory gown, household gloves and safety goggles. Prepare the working solution according to the dosing table. Run an ultrasonic cleaning cycle for 30 minutes. After cleaning, rinse the instruments thoroughly under running water of at least drinking-water quality for 10 seconds. Replace the working solution daily or when visibly contaminated. Adhere to the reprocessing recommendations of the instrument manufacturer. Not suitable for the cleaning of color-coded instruments or instruments made of natural bristles, nylon, or aluminum. A list of compatible and incompatible materials is available on our website. If you are unsure about the material composition of the medical device to be disinfected, contact the manufacturer before use. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**LIMITATIONS:** FOR PROFESSIONAL USE ONLY. Do not mix with other products. **WARNING:** DANGER H314 Causes corrosive to metals. H314 Causes severe skin burns and eye damage. P260 Do not breathe fume/mist/vapours/spray. P280 Wear protective gloves/protective clothing/eye protection/face protection. P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. P303 + P361 + P533 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P501 Dispose of contents/container in accordance with local and national regulations. Further information can be found in the Safety Data Sheet available on oroclean.com.

**STORAGE:** Store in a cool, dry and well-ventilated place. Keep away from direct sunlight, heat, and open flame. Storage temperature: < 30 °C. Dispose of the product after expiry date. **COMPOSITION:** 100 g contains 16 g potassium hydroxide, corrosion inhibitors, auxiliaries and water.

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