



Disinfectant wipes with a low alcohol content for medical devices.

- Aldehyde-, chlorine- and phenol-free
- Short exposure time
- High material compatibility
- Less than 50% alcohol content

actisolve® alcolite wipes are ready-to-use disinfectant wipes with a low alcohol content for the disinfection and cleaning of surfaces of non-invasive medical devices. The product can be safely used on sensitive medical devices such as treatment chairs with high-gloss surfaces. The disinfectant wipes are compatible with the most commonly used materials of medical devices as well as acrylic glass and artificial leather.

INSTRUCTIONS FOR USE

For your protection, wear disposable gloves and safety goggles. Place the refill pouch in the wipes container. Before use, cut open the refill pouch. Pull the first wipe at the heart of the roll through the star-shaped opening in the lid. Close the wipes box firmly after each use. Completely moisten and wipe the entire surface of the medical device. Leave disinfectant to act and adhere to contact time. Use on alcohol compatible materials only. Before first time use, test material compatibility by applying the disinfectant to a small inconspicuous area. A list of compatible and incompatible materials is available on our website. If you are uncertain as to the material composition of the medical device, contact the manufacturer before use. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

EFFECTIVE SPECTRUM / EXPOSURE TIMES

Bactericidal	60 sec
Antibiotic-resistant bacteria	60 sec
Yeasticidal	60 sec
Antibiotic-resistant yeasts	60 sec
Virucidal against enveloped viruses	15 sec
Bovine viral diarrhea virus (BVDV)	15 sec
Coronavirus	15 sec
Hepatitis B virus	15 sec
Hepatitis C virus	15 sec
Herpes simplex virus	15 sec
Human immunodeficiency virus (HIV)	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Limited spectrum virucidal	5 min
Adenovirus	5 min
Norovirus	5 min

Complies with EN 13624, EN 13727, EN 14476, and EN 16615.

COMPATIBILITY

Acetal polymer	Compatible
Acrylic glass	Compatible
Acrylonitrile butadiene styrene	Compatible
Artificial leather	Compatible
Chlorinated polyvinyl chloride	Compatible
Latex	Compatible
Polyamide	Compatible
Polybutylene terephthalate	Compatible
Polycarbonate	Compatible
Polyetherimide	Compatible
Polyethylene	Compatible
Polyethylene terephthalate	Compatible
Polyphenylene oxides	Compatible
Polyphenylsulfone	Compatible
Polypropylene	Compatible
Polystyrene	Compatible
Polysulfone	Compatible
Polyvinyl chloride	Compatible
Polyvinylidene fluoride	Compatible
Silicone	Compatible
Thermoplastic elastomer	Compatible
Thermoplastic polyurethane	Compatible
Ultra-high-molecular-weight polyethylene	Compatible

HANDLING

DANGER H228 Flammable solid. H319 Causes serious eye irritation. P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P233 Keep container tightly closed. P280 Wear protective gloves/protective clothing/eye protection/face protection. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P337 + P313 If eye irritation persists: Get medical advice/attention. P501 Dispose of contents/container in accordance to local and national regulations. For more information, please refer to Safety Data Sheet available on our website.

COMPOSITION

100 g moisturizing solution contains 7.5 g ethanol, 15 g isopropanol, 0.76 g quaternary ammonium compounds, auxiliaries and water.

PHYSICAL PROPERTIES

Physical state:	Liquid absorbed on fleece
Colour:	Colourless
Scent:	Of alcohol
pH-value:	Neutral
Shelf life:	24 months (unopened container)

DELIVERY FORMS

Saver, Refill pouch (14.5 x 20 cm), 250 pcs, classic scent
Soft, Refill pouch (19.3 x 21 cm), 120 pcs, classic scent

LEGAL MANUFACTURER



United Disinfectant Manufacturers AG
Dr. Grass-Strasse 12
9490 Vaduz
Principality of Liechtenstein
Telephone: +423 237 15 03
info@udm.li
www.udm.li